



INJECTION PROTOCOL

Allergenic extract prescription product for animal hyposensitization



	Dose	Day	Interval (Days)
Vial A Green	0.1 ml	1	0
	0.2 ml	3	2
	0.4 ml	5	2
	0.6 ml	7	2
	0.8 ml	9	2
	1.0 ml	11	2
	1.0 ml	13	2
	1.0 ml	15	2
	1.0 ml	20	5
	Vial B Blue	0.1 ml	25
0.2 ml		30	5
0.4 ml		36	6
0.6 ml		42	6
0.8 ml		48	6
1.0 ml		58	10
1.0 ml		68	10
1.0 ml		78	10
1.0 ml		88	10
Vial C Red		0.3 ml	98
	0.5 ml	112	14
	0.5 ml	133	21
	0.6 ml	163	30
	0.8 ml	193	30
	1.0 ml	223	30
	1.0 ml	255	30
CALL YOUR VET TO ORDER MAINTENANCE			
	1.0 ml	283	30

RESTRICTED TO USE BY OR UNDER THE SUPERVISION OF A VETERINARIAN

Contains phenol (0.4%) as a preservative

INDICATIONS AND USAGE

Keep treatment sets refrigerated at 2-7°C. Do not freeze. Injections are administered subcutaneously using sterile needle and syringes. **Contents should be brought to room temperature before injection.**

NEVER INJECTION INTRAVENOUSLY

Before injecting, withdraw the syringe plunger slightly. If blood appears in the syringe, reinsert in another site. The hyposensitization program may vary with each individual animal. The following injection schedule serves as a guideline. The actual dosage should always be determined by the veterinarian.

HANDLING ADVERSE REACTIONS, IF OBSERVED

Animals should be observed for at least 30 minutes following an injection. Possible systemic reactions to watch for are: increased itching, vomiting and/or diarrhea and hives. The possibility of an anaphylactic reaction is remote. However, if such a reaction occurs, it may be advisable to administer epinephrine. For an average size dog, injection 0.2 to 0.5ml of epinephrine 1:1000 w/v. Mild reactions may be controlled with antihistamines.

ADMINISTERING TWO TREATMENT SETS TO ONE ANIMAL

Using the same schedule, administer the recommended volume for each set preferably at different subcutaneous sites on the animal.

Please visit us on the web at: www.vetallergy.com
Print a personalized injection calendar for your patient.

Prescription products sales and support provided to **LICENSED VETERINARIANS ONLY**

Veterinarians call 800-553-1391 to order, and for TECHNICAL ASSISTANCE

SPECTRUM LABS

2143 E 5TH ST

TEMPE, AZ 85281

UNITED STATES VETERINARY LICENSE NUMBER 403